

TransZol Plant

Cat.No.: ET121

Storage: at room temperature for one year

Description

TransZol Plant is a ready-to-use reagent for the isolation of total RNA from polysaccharide-rich and/or polyphenol plant tissues, such as champignon, banana fruit, mango fruit, potato, carrot, sansevieria. It uses a modified CTAB method to lyse samples and phenol/chloroform to remove proteins and others impurities. It is also suitable for the isolation of total RNA from animal tissues like fat, connective tissues etc.

Highlights

- Superior lysis capability and higher RNA yield.
- The whole procedure can be completed in one hour.
- Pink solution for easy visualizing different phases.
- Unique dissolving solution for long-term RNA storage.

Kit Contents

Component	ET121-01
TP I buffer	100 ml
TP II buffer	100 ml
RNA Dissolving Solution	15 ml

Procedures

Reagents provided by customers: chloroform, isopropanol, 75 % ethanol (prepared with RNase-free water). All centrifugation steps are performed at 2-8°C.

- 1. Completely grind plant tissues in liquid nitrogen. Then add TP I solution to the homogenized tissues (add 1 ml of TPI for every 80-100 mg plant tissues). Mix thoroughly by pipetting several times and transfer all the lysate to a RNase-free microcentrifuge tube.
- 2. Centrifuge the lysate at 12,000×g for 5 minutes.
- 3. Transfer the supernatant (maybe a little cloudy) into two microcentrifuge tubes (about 400-500 μl each).
- 4. Add equal volume of TP II solution (pink) to each tube. Mix thoroughly by pipeting several times. Add chloroform (equal to 1/4 volume of the supernatant from step 3, about 100 -125 μl) to each tube, mix thoroughly by pipetting several times. Incubate at room temperature for 5 minutes.
- 5. Centrifuge the lysates at 12,000×g for 5 minutes. The lysates are separated into three layers: clear layer (colorless), interphase (colorless transparent oil, about 50 μl solution) and organic layer (pink). RNA is located in the clear layer.
- 6. Carefully transfer the colorless supernatant from two tubes to a new 1.5 ml RNase-free tube. (The clear layer and intermediate layer are difficult to distinguish. Recommend to leave about 50-100 µl colorless solution in the tube to avoid contamination).
- 7. Add equal volume of isopropanol to the transferred supernatant. Mix by inverting 4-6 times. Incubate at room temperature for 10 minutes.
- 8. Centrifuge at 12,000×g for 10 minutes. Discard the supernatant. RNA precipitate (white) can be seen on the bottom of the tube.
- 9. RNA pellet is dissolved in 30-40 µl of dissolving solution. For long-term storage, store the purified RNA at -70°C.

FOR RESEARCH USE ONLY